

---

# ECE 230L - Lab 3 Pre-lab

## ALICE Windows Installation


---

### Contents


|   |  |   |
|---|--|---|
| 1 | Windows 64 and 32 bit Installation Instructions  | 2 |
| 2 | Windows 32 bit Installation Further Instructions | 6 |
| 3 | Troubleshooting                                  | 8 |

# 1 Windows 64 and 32 bit Installation Instructions


1. Download the 64 Bit Libsmu software library installer package that contains the 64 Bit USB drivers for the ALM1000 [here](#). Download the link for the 32 bit libsmu library required to run ALICE1.3 for Windows [here](#).
2. Now we will install the USB drivers. Run the libsmu-1.0.2-setup-x64.exe program you just downloaded for 64 bit Windows and run the libsmu-1.0.2-setup-x86.exe program you just downloaded for 32 bit Windows. Windows will warn you that the software is from an unknown publisher. Say Yes anyway. You will first be prompted to select a Setup Language.




3. Accept the agreement below.




4. Next you will be asked to select a Destination Location. The default is fine. Click on Next.




5. Next you will be asked to select what you want to install. By default just the install WinUSB drivers will be selected. That is all you will need.




6. The next thing to pop up is just to confirm that you wish to install the software. Click on Install.




7. The Windows Driver install dialog box will now pop up. Click Next.



8. If the USB driver Install completes you will get the last pop up screen:




9. Click on Finish. The final pop-up should appear. Click on Finish.




## 2 Windows 32 bit Installation Further Instructions


1. Download the latest version of the Windows executable installer, alice-desktop-1.3-setup.exe, for ALICE 1.3 for the ALM1000 [here](#).
2. Next we will run the ALICE 1.3 Windows executable installer. Run the alice-desktop-1.3-setup.exe program you just down loaded. The first dialog to pop up asks you to select where you want to install the software. The default location is the best choice because the software writes certain files to the location where it runs so it needs to be somewhere that is not write protected. Click on Next.




3. The next pop-up asks you to select the Start Menu Folder. The default location is fine. Click on Next.



- The next pop-up just asks you to confirm that you want to install the software. Click on Next.



- When the installation completes this final dialog should appear. Click on Finish.




- After the software is installed a number of desktop icons should appear. You should now be able to plug-in your ADALM1000 to a USB port. Windows should recognize the new hardware and find the appropriate driver. You are now ready to run ALICE 1.3. Start ALICE by double clicking on the desktop Icon.




### 3 Troubleshooting

1. Most likely your ADALM1000 has firmware version 2.16 installed from the factory (unless you used Pixelpulse2 to update the firmware already). ALICE will run with this version of the firmware but a couple of advanced capabilities of the hardware will not be available. Starting ALICE 1.3 with 2.16 firmware will give you the screen below at start-up. You can press OK to continue.




2. This will happen each time ALICE starts until the new version 2.17 or higher of the firmware is flashed to the hardware. ALICE will then ask if you want to update the current firmware. Again you can skip this by saying No to use as is. This will happen each time ALICE starts until newer firmware is loaded. This firmware checking can be set to ignore by changing this line in the `alice_init.ini` file: `global IgnoreFirmwareCheck; IgnoreFirmwareCheck = 0` # change 0 to 1 to ignore firmware rev level.







- To update your board to the newest firmware click on the Options drop down menu. At the bottom of the list of options click on the Update Firmware button.




- The pop up dialog below will appear. Click on Yes.




- A file selection dialog will appear. You may need to navigate to where the software was installed C:\ALM Software\M1K. Select the m1000-2.17.bin file and click on Open.




- The program will pause for a few seconds while the board is being flashed. When it finishes the pop-up below should appear. Click on OK.




- You now need to unplug the board from the USB port to cycle the power. You also need to exit (close) ALICE. Click on OK and ALICE will close down.



- Plug the ADALM1000 back into the USB port. Restart ALICE by double clicking on the desktop Icon. ALICE should start normally now. You can check the current state of the board firmware by clicking on the About button under the File drop down menu.



9. If for some reason the flashing of the firmware fails to happen properly you can try again by clicking on the Update Firmware button again. If all else fails you will either need to install Pixelpulse2 and flash the firmware with that program or go back and fully install the libsmu library that includes the smu command line utility (selecting the Add application directory to path variable option). Using the Windows command prompt screen navigate to the directory where ALICE was installed C:\ALM Software\M1K and type smu. If the library was properly installed this should appear:



```
Command Prompt
C:\ALM Software\M1K>smu
smu: utility for managing M1K devices

-h, --help                print this help message and exit
--version                 show libsmu version
-l, --list-devices         list supported devices currently attached to the system
-p, --hotplug-devices      simple session device hotplug testing
-s, --stream-samples       stream samples to stdout from a single attached device
-d, --display-calibration  display calibration data from all attached devices
-r, --reset-calibration    reset calibration data to the defaults on all attached devices
-w, --write-calibration <cal file> write calibration data to a single attached device
-f, --flash <firmware image> flash firmware image to a single attached device

C:\ALM Software\M1K>smu -f m1000-2.17.bin
smu: successfully updated firmware
Please unplug and replug the device(s) to finish the process.

C:\ALM Software\M1K>
```

Using the `-f` command line option for smu should flash the new firmware onto the board in all cases. You may need to cycle the power to the board before it can be recognized and in be in the right mode to be flashed.